


Basic information

User Name: User1Date: 09-20-2023-17-05-52Total time: 0.58minRPM: 69.97

Height: --cmWeight: --KgGender: --Age: --yr

Motion Type Result

Knee Motion	Right	Left
Motion Type	Normal	Normal
Medial Percent	0.0%	0.0%
Normal Percent	100.0%	100.0%
Lateral Percent	0.0%	0.0%




Your right leg is Normal type cycling and left leg is Normal type cycling


Trajectory Result

Knee	Right	Left
Min Angle and Phase	-3.62 degree in Phase 4	-3.22 degree in Phase 4
Max Angle and Phase	1.38 degree in Phase 2	-0.59 degree in Phase 2

Right




Left




Your right leg motion range is about 5.0 degree and left leg is about 2.63 degree

Power Balance

Right



Left



Cycling Stable Result

Range 0-360	Right	Left
Stable	73.94%	68.86%
Average Std Angle	0.6 degree	0.65 degree

